

BRISTOL-PLUS™ | SURGICAL

PORTABLE LED HEADLIGHT SYSTEM



SEE BETTER
SEE FARTHER
SEE WHY

SINCE 1971
ENHANCING
YOUR VISION

Powered by
Long Island
Technology
Group


BRISTOL-PLUS™ SURGICAL

MEDIUM-INTENSITY, PORTABLE LED SURGICAL HEADLIGHT WITH ADJUSTABLE SPOT, POWERED BY A 3-CELL BATTERY

- Illumination Range: ~40,000 LUX at 16" (40 cm) to ~132,000 LUX at 10" (25 cm)
- Adjustable Spot Diameter: 1"-3" (2.5 - 7.6 cm) at a working distance of 16" (40 cm)
- Color Rendering Index: 78 CRI
- Color Temperature: 6,500 K
- Intensity Settings: 5 Levels
- 3-Cell Battery
- Battery Duration at Full Intensity: 10 hours
- Battery Expected Life Span: > 300 Charges
- LED Expected Life Span: > 10,000 hours
- Optic Weight: 1.2 oz (34 g)
- Battery Weight: 6 oz (170 g)



High Brightness at all Working Distances (WD)



BRISTOL-PLUS PORTABLE LED HEADLIGHT SYSTEM INCLUDES:

- Bristol-Plus Surgical Headlight
- Choice of: Fiji or Key West Headband, or Optic Clip
- Choice of One or Two 3-Cell Battery
- Charging Cradle, Power Supply
- Holster, Linkage Kit, and Carrying Case

BRISTOL-PLUS PORTABLE LED HEADLIGHT PRODUCT CODES:

BFW 9040P	Bristol-Plus Headlight System with Fiji Headband with One 3-Cell Battery
BFW 9040PB	Bristol-Plus Headlight System with Fiji Headband and Two 3-Cell Batteries
BFW 9042P	Bristol-Plus Headlight System with Optic Clip with One 3-Cell Battery
BFW 9042PB	Bristol-Plus Headlight System with Optic Clip and Two 3-Cell Batteries
BFW 9044P	Bristol-Plus Headlight System with Key West Headband with One 3-Cell Battery
BFW 9044PB	Bristol-Plus Headlight System with Key West and Two 3-Cell Batteries
BFW 1552C	Bristol-Plus Optic Clip
BFW 9022ACC	Bristol-Plus Charger
BFW 9035	Bristol-Plus 3-Cell Battery
BFW 9036	Bristol-Plus Holster

EXCELLENT FOR USE IN THE FOLLOWING DEPARTMENTS

HOSPITAL

- ENT
- GENERAL SURGERY
- CARDIOSURGERY
- NEUROSURGERY
- ORTHOPEDIC SURGERY
- PEDIATRIC SURGERY
- SPINE SURGERY
- THORACIC SURGERY
- TRANSPLANT SURGERY
- VASCULAR SURGERY

SURGERY CENTER

- ORTHOPEDIC SURGERY
- PLASTIC SURGERY

OFFICE

- DENTAL EXAM
- ORAL/MAXILLOFACIAL



BFWinc.com

800-717-4673 | +1-502-899-1808

445 BAXTER AVENUE | SUITE 175

LOUISVILLE, KY 40204 USA



BRISTOL-PLUS TESTED AND CERTIFIED TO:
IEC 60601-1 Edition 3, IEC 60601-1-2 Edition 4
ANSI/AAMI ES60601-1 (2005 + C1:09 + A2:10),
CAN/CSA-C22.2 NO 60601-1- (2008), CE

ENGINEERED AND ASSEMBLED IN THE U.S.A. © 2020 BFW, Inc.
ALL TECHNICAL SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE